Tubes Radiogènes à Anode Tournante  
 Drehanoden - Röntgenröhre  
 Tubos de Rayos-X con Ánodo Giratorio

## GS-3576HS/MX165CT Housing

Note: Document originally drafted in the English language.

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

Manufactured by Varian Medical Systems  
Fabrique par Varian Medical Systems  
Hergestellt von Varian Medical Systems  
Fabricado por Varian Medical Systems

Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.

Caractéristiques d'Émission du Filament CEI 60613  
Kathoden - Emissionskennlinien IEC 60613  
Características de Emisión del Catodo IEC 60613

Le Gaine MX165CT

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Encaje de MX165CT

Maximum Peak Voltage .....	140 kV
Anode to Ground .....	70 kV
Cathode to Ground .....	70 kV
Maximum X-ray Tube Assembly Heat Content .....	2.6 MJ (3.7 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat) .....	3.7 kW (5.0 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly Permanent filtration .....	1.2 mm Al IEC 60522
Loading Factors for Leakage Radiation .....	140 kV, 26 mA
Ambient Air Temperature Limits for Operation .....	5°C to 35°C
Temperature Limits for Storage and Transport .....	-34°C to +60°C
Humidity .....	10% to 90%
Atmospheric Pressure Range .....	70 kPa to 106 kPa
Weight - Tube and Housing .....	75.6 kg (167 lbs)
IEC Classification .....	Class 1
Safety Devices - Thermal Switch, Normally Closed Contact .....	Opening at 74°C to 79°C
Pressure Switch Normally Closed Contact .....	Opening at 5 PSIG ±3 PSIG

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
**GS-3576HS**  
**TUBE HOUSING ASSEMBLY HEATING AND COOLING**

**Note:**

- Heat inputs into housing include tube power, filament power, and stator power.
- Heating curves based on no restrictions of natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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


<b>Stator:</b>		
<b>Stator Coil Resistance:</b>		
Black to White	23 Ohms ±15%	
Green to White	23 Ohms ±15%	
Black to Green	46 Ohms ±15%	
<b>Starter Voltage:</b>		
	<b>Start</b>	<b>Run</b>
50/60 Hz	385 VAC	75 VAC
150/180 Hz	450 VAC	85 VAC
<b>Time to Full Speed:</b>		
50/60 Hz	0 - 2700 RPM	10 Sec.
<b>X-Ray Tube Assembly:</b>		
GS-3576HS/MX165CT IEC 60601-2-28		

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
## ANODE HEATING AND COOLING CURVES

**Note:**

1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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