Tubes Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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Manufactured by Varian Medical Systems
Fabrique par Varian Medical Systems
Hergestellt von Varian Medical Systems
Fabricado por Varian Medical Systems

Specifications subject to change without notice.
Spécifications susceptibles d'être modifiées sans préavis.
Technische Daten ohne Gewähr.
Especificaciones sujetas a cambio sin previo aviso.

Tableaux des Caractéristiques Nominales de Balayage Volumétrique/Hélicoïdal CEI 60613
Volumen-/Spiralbelichtungs-Leistungdiagramme IEC 60613
Volumétrico/Clasificación Grafica del Escán/Helicoideo IEC 60613

3 Ø 180 Hz ■

0.8 x 0.4 Focal Spot 7 Degrees
0,8 x 0,4 Dimension Focale 7 Degrés
0.8 x 0.4 Brennpunkt 7 Grad
0.8 x 0.4 De Marcas Focales 7 Grados

| Volume Scan Time (Seconds) | MAXIMUM ALLOWED TUBE CURRENT (mA) | | | | | | | | |
|----------------------------------|--|--------|--------|----------------------|--------|--------|----------------------|--------|--------|
| | AS A FUNCTION OF THE FOLLOWING STARTING HEAT STORAGE AND TUBE VOLTAGES | | | | | | | | |
| | Starting H.S. = 40 % | | | Starting H.S. = 55 % | | | Starting H.S. = 70 % | | |
| | 100 kV | 120 kV | 130 kV | 100 kV | 120 kV | 130 kV | 100 kV | 120 kV | 130 kV |
| 2 | 190 | 160 | 140 | 190 | 160 | 140 | 190 | 160 | 140 |
| 3 | 190 | 160 | 140 | 190 | 160 | 140 | 190 | 160 | 140 |
| 4 | 190 | 160 | 140 | 190 | 160 | 140 | 190 | 160 | 140 |
| 5 | 190 | 160 | 140 | 190 | 160 | 140 | 190 | 160 | 140 |

Note:
1. Limits are based on maximum track rating except for the following codes:
a - Limited by available heat storage.
b - Limited by window heating.
c - Limited by filament emission.
2. H.S. = Heat Storage
kV = Tube Voltage

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Note:
Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

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
Le Gaine B-172H/MX115CT

Das B-172H/MX115CT Gehäuse

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| | |
|---|-------------------------|
| Maximum Peak Voltage | 140 kV |
| Anode to Ground | 70 kV |
| Cathode to Ground | 70 kV |
| Maximum X-ray Tube Assembly Heat Content | 2.2 MJ (3.0 MHU) |
| Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat) | 4.0 kW (340 kHU/min) |
| X-Ray Tube Assembly - IEC 60522 | |
| Permanent filtration | 1.2 mm Al at 140kV |
| Collimator | 1.5 mm Al at 140kV |
| Total Filtration | 2.7 mm Al at 140kV |
| Loading Factors for Leakage Radiation | 140 kV, 21 mA |
| Federal Standard High Voltage Cable | 72 |
| Ambient Air Temperature Limits for Operation | 5°C to 40°C |
| Temperature Limits for Storage and Transport | -20°C to +75°C |
| Humidity | +10% to +90% |
| Atmospheric Pressure Range | 70 kPa to 106 kPa |
| Weight: Housing and Heat Exchanger | 52.6 kg (116 lbs) |
| IEC Classification | Class I |
| Safety Devices - Thermal Switch | |
| Normally Closed Contact | Opening at 74°C to 79°C |
| Pressure Switch | |
| Normally Closed contact | Opens at 5PSIG ± 3PSIG |
| Power Supply | DC |

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

Note:

- Heat inputs into housing include tube power, filament power, and stator power.
- Heating curves based on no restrictions of natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

| | | |
|---|----------------|------------|
| Stator Type: GE or SARC | | |
| Stator Coil Resistance: | | |
| Black to White | 23.5 Ohms ±15% | |
| Green to White | 23.5 Ohms ±15% | |
| Starter Voltage: | Start | Run |
| 150/180 Hz | 385 VAC | 85 VAC |
| Time to Full Speed: | | |
| 150/180 Hz | 0 - 7500 RPM | 4.4 Sec. |
| X-Ray Tube Assembly: | | |
| GS-1562/B-172H | | |
| GS-1562/MX115CT IEC 60601-2-28 | | |
| (Test based on "Maximum Continuous Anode Input Rate") | | |

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Note:
1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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